

BOA, IV Constricting Band



Item No: CL-0018

Manufacturer: North American

Rescue Products

Certificates: CE, Latex Free

Red Tip Technology

19 to 45cm Lenghts available:

25 to 68cm

Country of origin: USA

NSN no: 6515-01-537-2611

Tariff: 90189084

Weight: 42q Increase the proficiency in gaining direct IV access with the BOA® Constricting Band, a now-critical component of successful IV cannulation.

Not since constricting bands were first conceived has the design been so fundamentally altered as it has been with the BOA®.

For improved central venous access every time, the BOA® uses a rolling action to transfuse venous blood to the distal veins while still allowing arterial flow. (For maximum effect, the BOA® is placed as high up on the extremity as allowable.) This technique distends the veins to make IV access quicker and more effective even in the most challenging pa-

The BOA®'s patented design features a Quick Click connector and two-finger Quick Release for added efficiency. No longer do healthcare professionals need to tie (and often re-tie) a constricting band. The BOA® decreases "time on task" for treating combat casualties, civilian trauma patients, or medical patients with poor vascular access.

7550 Sørvad Denmark